

# **ALPHA** Unpressurised vessel series for small amounts of nitrogen.



Safe and lightweight aluminum storage and transport vessel for liquid nitrogen.



00:A

100

9

-250

- 200

- 150

50:1 ml 50

40

- 30

20-

## All Advantages at a Glance

- Small, Compact and Sturdy
- Available in three sizes with a capacity of approx. 6, 10 & 25 litres
- Comprehensive accessories



08

60

www.cryotherm.de | www.cryotherminc.com

## ALPHA

## **Technical Data**

Storage Capacity ALPHA		5	10	25
Net capacity	[]	6	10	25
Daily evaporation	[l/Day]*	0.12	0.12	0.12
Static holding time	[Days]*	50	83	208
Weight, empty	[kg]	5	6.2	10.5
Weight, full (filled with LIN)	[kg]	9.7	13.9	31.6
Total height	[mm]	468	528	700
Neck tube diameter	[mm]	50	50	50
Inner diameter	[mm]	231	231	296
Outer diameter	[mm]	300	300	396

لللالالال

\*) The daily evaporation rate and the static holding time were determined at 20 degrees Celsius, 1013 hectopascal, with the vessel standing still, lid closed and without accessories. These are nominal values that may change depending on the use and age of the container.

### Suitable accessories



 Cryo-Gloves

 Length: 400 mm

 Item No.

 Size S:
 78400843

 Size M:
 78400180

 Size L:
 78400887

 Size XL:
 78400856



Protection shield Item No. 77031452



Roller base for ALPHA 25 with five swivel castors, three of them lockable. Item No. 78217449

Size XL: 784008



We supply various equipment variants and a wide range of accessories for your ALPHA.

<b>ALPHA Lines</b>		
🗴 Item No.		Quantity
78224172	ALPHA 5, pressureless transport/storage vessel for liquid nitrogen	
78224173	ALPHA 10, pressureless transport/storage vessel for liquid nitrogen	
78224174	ALPHA 25, pressureless transport/storage vessel for liquid nitrogen	
78217449	Roller base for ALPHA 25, with five swivel castors, three of them lockable	
Other container	sizes on request.	

sales@cryotherm.de

**\$** +49 (0)2741 - 95 85 0

#### Cryotherm GmbH & Co. KG

Euteneuen 4 57548 Kirchen/Sieg, Germany Telephone: ++49 27 41 / 95 85-0 Fax: ++49 27 41 / 69 00 E-mail: info@cryotherm.de Version 2032A



#### Cryotherm Inc.

5126 South Royal Atlanta Drive Tucker, GA 30084, USA Telephone: +1 770-274-2410 Fax: +1 770-934-3384 E-mail: info@cryotherminc.com







Our BIOSAFE® ALPHA serie is available with a vessel capacity of 6 to 25 litres.

Further infor-

mation in the Lifescience catalogue

CNE

Life Science